



## O2 CONCENTRATOR

### QUALITY

- ✓ Triple air filtration including carbon and bacterial filtration.
- ✓ Multi parameter display, including percentage of concentration.
- ✓ Easy filter access.
- ✓ Consistent supply of high concentration supplemental oxygen.

### ECONOMICAL

- ✓ Has the ability to pay for itself in less than a year.
- ✓ Low total cost of ownership.
- ✓ Eliminates the requirement of delivered oxygen.

### SAFE

- ✓ Power failure, no flow, high and low pressure alarms.

### CONVENIENT

- ✓ Smooth quiet operation.
- ✓ Casters for portability.
- ✓ Programable digital run timer.
- ✓ Easily adaptable to any anesthesia system (O2 DISS MALE and Hose Barb fitting included)

### SPECIFICATIONS

Weight	5 lpm: 60lbs OR 10 lpm: 65lbs
Height	25 in
Depth	12.5 in
Width	15 in
Flow Meter Capacity	0-5 lpm OR 0-10 lpm
Sound	40 - 45 dBa
Line Pressure	10 psi

1-877-637-3625

[info@somniscientific.com](mailto:info@somniscientific.com)

[www.somniscientific.com](http://www.somniscientific.com)

### INHALATION ANESTHESIA EQUIPMENT AND SERVICE

**SOMNI Scientific** is centered on the animal welfare and research community with a focus on clinical accuracy, clinician/technician safety, economic performance and intuitive functionality.

