



REBREATHING SYSTEM WORK STATION VETERINARY REBREATHING ANESTHESIA SYSTEM




QUALITY:

- H-frame provides exceptional support and stability
- 3" wheels for ease of transport
- High quality, durable, limited maintenance design
- Silicone tubing and gaskets throughout

▶ COMPATIBLE ACCESSORY OPTIONS

Dual Vaporizer
Bain Block Manifold
Dual Gas supply
Dual E-cylinder
O2 Concentrator mount
5 year warranty



SAFETY:

- Pressure Manometer center mount for immediate visibility
- Popoff valve releases excess pressure at 2cmH20 while maintaining passive volume in the rebreathing bag
- Toggle switch: Nonrebreathing circuit connection stays in place, no tubing to connect/disconnect



CONVENIENT:

- Lightweight, small footprint, easy assembly
- Easily Configured shelves for additional workspace or accessories
- Powder coated and anodized frame for easy cleaning



1-877-637-3625

info@somnivets.com

www.somnivets.com

Proudly machined, assembled, and calibrated in the USA