

SOMNI
SCIENTIFIC



SOMNI CATCH
STEREOTAXIC WASTE GAS
COLLECTION

SOMNI Catch

(Stereotaxic Waste Gas Collection)

The SOMNI Catch is a “V-Channel” Waste Anesthetic Gas scavenging device developed for use with most models of stereotaxic surgery platforms. The SOMNI Catch effectively removes the potentially harmful waste gases that can leak into the environment.



QUALITY:

- Machined from anodized aluminum for durability and chemical resistance.
- Large bore vacuum port for maximum airflow “V” shape design to channel waste gas more efficiently.


► **SPECIFICATIONS**

WEIGHT: 4.7 oz. / 146.19 g

HEIGHT: 1 in / 2.54 cm

DEPTH: 2 in / 5.08 cm

WIDTH: 2 in / 5.08 cm




SAFETY:

- Prevents unnecessary exposure to personnel and environment (NIOSH 2ppm).
- Extensively tested, effective design.



CONVENIENT:

- Ability to connect to SOMNI EPS-3 or in-house active scavenging system.
- Configurable to accommodate left or right side work station.



1-877-637-3625

info@somniscientific.com

www.somniscientific.com

Proudly machined, assembled, and calibrated in the USA