



# SOMNI BAIN MANIFOLD



## SAFETY:

- Pressure Manometer center mount for immediate visibility
- Popoff valve releases excess pressure at 2cm H<sub>2</sub>O while maintaining passive volume in the rebreathing bag
- Monitoring pressure during manual ventilation is critical in preventing barotrauma or inadequate ventilation, especially in small patients

Standard anesthesia systems do not incorporate this feature

## ► SPECIFICATIONS


WEIGHT: 2 lbs.

HEIGHT: 8 in.

DEPTH: 7.5 in.

WIDTH: 3 in.

SOMNI Bain Nonrebreathing Assembly shown with RM-17003 Small Animal Rebreathing Anesthesia Machine and VS-6001 SOMNI 3 Sevoflurane Vaporizer Funnel Fill 250 ML Capacity



## CONVENIENT:

- Simple press and twist operation to anchor.
- Freely rotates on the base.
- Easy and quick to move when flipping the patient.



## QUALITY:

- Manufactured from high quality, durable materials
- Adjustable manometer
- Space conscious
- Easily adapted to most breathing systems



1-877-637-3625

info@somnivets.com

www.somnivets.com

Proudly machined, assembled, and calibrated in the USA