

SOMNI
SCIENTIFIC



SOMNI BAIN MANIFOLD



SAFETY:

- Pressure Manometer center mount for immediate visibility.
- Popoff valve releases excess pressure at 20cm H₂O while maintaining passive volume in the rebreathing bag.
- Monitoring pressure during manual ventilation is critical in preventing barotrauma or inadequate ventilation, especially in small patients.

- Standard anesthesia systems do not incorporate this feature.


► SPECIFICATIONS

WEIGHT: 2 lbs.

HEIGHT: 8 in.

DEPTH: 7.5 in.

WIDTH: 3 in.




CONVENIENT:

- Simple press and twist operation to anchor.
- Freely rotates on the base.
- Easy and quick to move when flipping the patient.



QUALITY:

- Manufactured from high quality, durable materials.
- Adjustable manometer
- Space conscious
- Easily adapted to most breathing systems.



1-877-637-3625



info@somniscientific.com



www.somniscientific.com

Proudly machined, assembled, and calibrated in the USA