

SOMNI
SCIENTIFIC



O2 CONCENTRATOR



QUALITY:

- Triple air filtration including carbon and bacterial filtration.
- Multi parameter display, including percentage of concentration.
- Easy filter access.
- Consistent supply of high concentration supplemental oxygen.

SPECIFICATIONS

WEIGHT: 5 lpm: 35 lbs OR 10 lpm: 38 lbs

Height: 25 in / 635 mm

Depth: 12.5 in / 317.5 mm

Width: 15 in / 381 mm

FLOW METER CAPACITY: 0-5 lpm OR 0-10 lpm

SOUND: 40 - 45 dBa

LINE PRESSURE: 10 psi

POWER OPTIONS: Domestic 110V, International 220V



ECONOMICAL:

- Has the ability to pay for itself in less than a year.
- Low total cost of ownership.
- Eliminates the requirement of delivered oxygen.




SAFE:

- Power failure, no flow, high and low pressure alarms.



CONVENIENT:

- Smooth quiet operation.
- Casters for portability.
- Programmable digital run timer.



1-877-637-3625

info@somniscientific.com

www.somniscientific.com

Proudly machined, assembled, and calibrated in the USA