

# SOMNI BAIN NON-REBREATHING ASSEMBLY

## QUALITY

- ✓ Manufactured from high quality, durable materials
- ✓ Adjustable manometer
- ✓ Space conscious
- ✓ Easily adapted to most rebreathing anesthesia systems

## SAFETY

- ✓ Pressure Manometer center mount for immediate visibility
- ✓ Popoff valve releases excess pressure at 2cmH<sub>2</sub>O while maintaining passive volume in the rebreathing bag
- ✓ Monitoring pressure during manual ventilation is critical in preventing barotrauma or inadequate ventilation, especially in small patients  
Standard anesthesia systems do not incorporate this feature

## CONVENIENT

- ✓ Included CPRAM allows enough length for imaging and specialty procedures
- ✓ Tube within tube design allows warming & humidification of gases
- ✓ Waste gas scavenging port

## SPECIFICATIONS

Weight	3 lbs.
Height	8 in
Depth	7.5 in
Width	3 in



SOMNI Bain Nonrebreathing Assembly shown with RM-17003 Small Animal Rebreathing Anesthesia Machine and VS-6001 SOMNI 3 Sevoflurane Vaporizer Funnel Fill 250 ML Capacity

**SOMNI PROVIDES UNPARALLELED CUSTOMER SERVICE, CLINICAL AND TECHNICAL SUPPORT.**

For more information or pricing, please contact our sales team.

☎ 1-877-637-3625    ✉ info@somniscientific.com

🌐 www.somniscientific.com



*Proudly machined, assembled  
and calibrated in the USA*