

SOMNI
SCIENTIFIC



SOMNI REBREATHING ANESTHESIA SYSTEM WALL MOUNT STYLE



QUALITY:

- Low Profile mount to save space in OR.
- Multi piece design makes for easy mounting.
- High quality, durable, limited maintenance design.
- Silicone tubing and gaskets throughout.
- Designed to switch from a rebreathing circuit to a non-rebreathing circuit with a simple toggle switch.

COMPATIBLE ACCESSORY OPTIONS

Bain Block Manifold

Dual Gas supply

5 year warranty



SAFETY:

- Pressure Manometer center mount for immediate visibility.
- Popoff valve releases excess pressure at 20cm H2O while maintaining passive volume in the rebreathing bag.
- The oxygen toggle switch works as a built-in safety system, that keeps the pressure at 20PSI. This is designed with patient and user safety in mind.



CONVENIENT:

- Save space by wall mounting!
- 90 degree gas fittings on back for easy O2 connections.
- Anodized frame for easy cleaning.



1-877-637-3625

info@somniscientific.com

www.somniscientific.com

Proudly machined, assembled, and calibrated in the USA