

AMD-3

The Somni AMD-3 accurate multiple delivery is a flowmeter controlled diverter system that is adaptable to any standard Inhalant Anesthesia System. It can be added to Somni's TT-2 or any other non-pressurized system. The AMD-3 gives the user the ability to delivery accurate flowrates to 3 stations simultaneously. No more guessing if the flow is being delivered equally to multiple devices at once. With its included pressure relief valve the user does not have to worry about pressure building up in the vaporizer. 3 fresh gas common outlets supplied with color coded tubing makes this system very user friendly.

FEATURES

- ✓ 3 Flowmeter Controlled Stations with Toggle Switches
- ✓ Color Coded Tubing to each station
- ✓ Common Gas outlets for each station
- ✓ Small footprint
- ✓ Adjustable Flowrates
- ✓ Stacking design saves valuable workbench space

INCLUDES

- ✓ Pressure Relief Valve
- ✓ 3 each, 6' color coded tubing with adapters

SPECIFICATIONS

Weight	
Height	7 in / mm
Depth	8 in / mm
Width	12 in / mm



SOMNI PROVIDES UNPARALLELED CUSTOMER SERVICE, CLINICAL AND TECHNICAL SUPPORT.

For more information or pricing, please contact our sales team.

☎ 1-877-637-3625 ✉ info@somniscientific.com

🌐 www.somniscientific.com



*Proudly machined, assembled
and calibrated in the USA*